


# DNA ligation & transformation


TA. Oh Hyeokjin  
(basic science building 515)

# What we will do?

- PCR Amplification of DNA of interest (DOI)
- Purification of amplified POI
- Restriction of amplified POI and plasmid
- Purification of restricted DNAs
- Ligation of POI and plasmid
- Transformation of ligated DNA
- DNA sequencing
- Overexpression
- purification




# DNA ligation



Using DNA Ligase to Splice Together Sticky-Ended DNA Fragments

# DNA ligation




# DNA ligation

## Materials


DDW, T4 ligase ligation buffer, T4 ligase, restricted DNAs (target of interest and plasmid)

1. Thaw digested (restricted) DNA of interest (DOI) and plasmid and T4 ligase ligation buffer.
2. Prepare two new 1.7 ml micro-tubes.
3. Data from last week, calculate relative amount of DOI to plasmid.
4. Add DOI and plasmid into a micro-tube → insert ligation sample (Ratio of DOI and plasmid should be 3:1 ~ 5:1)  
Add DOI only into a micro-tube and add DDW instead of plasmid → self ligation sample
5. Add 1 ul of T4 ligase ligation buffer.
6. Add DDW up to 9 ul.
7. Add 1 ul of T4 ligase (total volume should be 10 ul).
8. Close lid of micro-tube and centrifuge, vortex and centrifuge.
9. Incubate for 2h at 37°C.

# Transformation




# Transformation




created by North Carolina School of Science and Mathematics

copyright © 2013 North Carolina Department of Public Instruction

 Except where otherwise noted, this work is licensed under <http://creativecommons.org/licenses/by-nc-sa/3.0/>

## Heat-shock transformation

# Transformation



electroporation transformation



# Heat-shock Transformation

## Materials

competent cell (DH5a), ligated DNA, 42°C heat-block, LB

1. Thaw competent cell **in ice**.
2. Prepare two new 1.7 ml micro-tubes.
3. Add 20 ul of competent cell into micro-tubes.
4. Add 2 ul of insert ligation sample and self ligation sample to each tube. (Generally, the ratio between cell and gene is 10 : 1.)
5. Tapping and incubate for 10 min **in ice**.
6. Move the micro-tube to 42°C heat-block and incubate for 30 sec.
7. After incubation, move the micro-tube **to ice** and incubate for 1 min.
8. Add 800 ul of LB medium.
9. Incubate for 30 min at 37°C.
10. Prepare LB-agar plate that contains appropriate antibiotics.
11. Spread the cell+gene solution in the LB-agar plate with spreader.
12. Incubate for 12~18h at 37°C.