

# Synthesis of $\text{Cr}(\text{acac})_3$ and $\text{Co}(\text{acac})_3$

# Process of Inorganic Chemistry Experiment

**Synthesis**

```
graph TD; A[Synthesis] --> B[Purification]; B --> C[Characterization];
```





**Purification**

- **Recrystallization**
- Column chromatography

**Characterization**


- **NMR**, IR, UV, Single Crystal X-ray diffraction
- Elemental Analysis, Mass Spectrometry

# Synthesis of $M(\text{acac})_3$ ( $M = \text{Cr}$ and $\text{Co}$ )




# Chelate Effect


→ Chelating ligands have a higher affinity to binding to the metal compared to the monodentate ligands.




$$\Delta G = -RT \ln K = \Delta H - T\Delta S$$




Monodentate



Bidentate



Tridentate



Tetradentate

# Procedure 1. Synthesis

## <Cr(acac)<sub>3</sub>>

- 50 mL of **distilled water** & 2.66 g of **chromium(III) chloride hexahydrate** in a 100 mL Erlenmeyer flask.
- Add 10 g of **urea** & 5.9 mL of **acetylacetone**.
- Heat the mixture for 1.5 hr.

## <Co(acac)<sub>3</sub>>

- 1.25 g of **cobalt(II) carbonate** and 10 mL of **acetylacetone** in a 50 mL Erlenmeyer flask.
- Heat the mixture to 100 °C.
- Transfer the flask from the water bath, and add 2.5 mL of 30 % **H<sub>2</sub>O<sub>2</sub> DROPWISE. (Caution: Vigorous reaction)**
- Reheat, then add more 2.5 mL of 30 % **H<sub>2</sub>O<sub>2</sub> DROPWISE.**
- Heat the mixture to 100 °C for 1 hr.

## •Tips!!

- **First heat the water bath and then weigh the compounds**
- **Put two flasks in the same boiling water bath.**
- **Prepare a cold ethanol for the purification**



# Procedure 1. Purification

## <Cr(acac)<sub>3</sub>>

- Cool down the reaction mixture to room temperature.
- Filter the mixture and wash the violet solid with distilled water.
- Dry and measure the mass.

## <Co(acac)<sub>3</sub>>

- Cool down the reaction mixture in an ice-salt bath for 30 minutes.
- Filter the mixture and wash the green solid with distilled water and small amount of cold ethanol.
- Dry and measure the mass.




# Procedure 2. Purification-Recrystallization

| Solvent                                   | Boiling point | Dielectric constant | Density           |
|---|---------------|---------------------|-------------------|
| <b>Non-Polar Solvents</b>                 |               |                     |                   |
| <a href="#">Hexane</a>                    | 69 °C         | 2.0                 | 0.655 g/ml        |
| <a href="#">Benzene</a>                   | 80 °C         | 2.3                 | 0.879 g/ml        |
| <a href="#">Toluene</a>                   | 111 °C        | 2.4                 | 0.867 g/ml        |
| <a href="#">Diethyl ether</a>             | 35 °C         | 4.3                 | 0.713 g/ml        |
| <a href="#">Chloroform</a>                | 61 °C         | 4.8                 | <b>1.498 g/ml</b> |
| <a href="#">Ethyl acetate</a>             | 77 °C         | 6.0                 | 0.894 g/ml        |
| <b>Polar Aprotic Solvents</b>             |               |                     |                   |
| <a href="#">1,4-Dioxane</a>               | 101 °C        | 2.3                 | <b>1.033 g/ml</b> |
| <a href="#">Tetrahydrofuran (THF)</a>     | 66 °C         | 7.5                 | 0.886 g/ml        |
| <a href="#">Dichloromethane (DCM)</a>     | 40 °C         | 9.1                 | <b>1.326 g/ml</b> |
| <a href="#">Acetone</a>                   | 56 °C         | 21                  | 0.786 g/ml        |
| <a href="#">Acetonitrile (MeCN)</a>       | 82 °C         | 37                  | 0.786 g/ml        |
| <a href="#">Dimethylformamide (DMF)</a>   | 153 °C        | 38                  | 0.944 g/ml        |
| <a href="#">Dimethyl sulfoxide (DMSO)</a> | 189 °C        | 47                  | <b>1.092 g/ml</b> |
| <b>Polar Protic Solvents</b>              |               |                     |                   |
| <a href="#">Acetic acid</a>               | 118 °C        | 6.2                 | <b>1.049 g/ml</b> |
| <a href="#">n-Butanol</a>                 | 118 °C        | 18                  | 0.810 g/ml        |
| <a href="#">Isopropanol (IPA)</a>         | 82 °C         | 18                  | 0.785 g/ml        |
| <a href="#">n-Propanol</a>                | 97 °C         | 20                  | 0.803 g/ml        |
| <a href="#">Ethanol</a>                   | 79 °C         | 24                  | 0.789 g/ml        |
| <a href="#">Methanol</a>                  | 65 °C         | 33                  | 0.791 g/ml        |
| <a href="#">Formic acid</a>               | 100 °C        | 58                  | <b>1.21 g/ml</b>  |
| <a href="#">Water</a>                     | 100 °C        | 80                  | 1.000 g/ml        |

## Procedure 2. Purification-Recrystallization

- Take ~ 0.2 g of each sample.
- Dissolve in warm **toluene**. (small amount) to form saturated solution.
- Add few drops of **petroleum ether slowly**. → crystals will form.
- Add more **petroleum ether**.
- Filter and wash with **petroleum ether**.
- Dry and measure mass .

<Cr(acac)<sub>3</sub>>



<Co(acac)<sub>3</sub>>





## Procedure 3. Characterization-UV

- The energy in the UV and visible light region is used to excite species to higher electronic energy levels.
- The absorbance, **A**, is defined as

$$\mathbf{A} = \log(\mathbf{I}_0/\mathbf{I}), \text{ where } \mathbf{I}_0 = \text{incident intensity} \\ \text{and } \mathbf{I} = \text{measured intensity}$$

- Beer-Lambert Law

$$\mathbf{A} = \boldsymbol{\varepsilon} \mathbf{b} \mathbf{c}, \text{ where } \boldsymbol{\varepsilon} = \text{extinction coefficient (L}\cdot\text{mol}^{-1}\cdot\text{cm}^{-1}), \\ \mathbf{b} = \text{optical path length (cm),} \\ \text{and } \mathbf{c} = \text{molar concentration (mol}\cdot\text{L}^{-1})$$

$\boldsymbol{\varepsilon} = 10^3 \sim 10^5 \text{ dm}^3\cdot\text{mol}^{-1}\cdot\text{cm}^{-1}$  for fully allowed transition.

$\boldsymbol{\varepsilon} = 10^0 \sim 10^3 \text{ dm}^3\cdot\text{mol}^{-1}\cdot\text{cm}^{-1}$  for orbitally forbidden transition.

$\boldsymbol{\varepsilon} = 10^{-5} \sim 10^0 \text{ dm}^3\cdot\text{mol}^{-1}\cdot\text{cm}^{-1}$  for spin forbidden transition.

# Procedure 3. Characterization-UV

## ***Selection Rule***

### 1. Laporte Selection rule

→ For centrosymmetric molecules, electronic transitions that conserve parity are **forbidden**.

(g → g or u → u)

### 2. Spin Selection rule

→ Transitions with same spin state is **allowed**.


ex) Metal d-d transition

→ Spin allowed, Laporte forbidden


# Procedure 3. Characterization-NMR

One is  $\text{Co}(\text{acac})_3$  and another is  $\text{Cr}(\text{acac})_3$ . Which is which?

(A)




(B)



# Procedure 3. Characterization-NMR

This is acetylacetone NMR data. Explain this peak splitting.



# Procedure 3. Characterization-NMR

d-orbital splitting of the octahedral complex

