

Exp.5


# Polymerase Chain Reaction (PCR)

# Biological Engineering and Cloning

**Biological engineering, or bioengineering/bio-engineering,** is the application of principles of biology and the tools of engineering to create usable, tangible, economically viable products.


**Cloning** is the process of producing genetically identical individuals of an organism either naturally or artificially. In nature, many organisms produce clones through asexual reproduction.

**Cloning in biotechnology refers to the process of creating clones of organisms or copies of cells or DNA fragments** (molecular cloning).



# Recombinant DNA

Genetically engineered DNA prepared by transplanting or splicing genes from one species into the cells of a host organism of a different species. Such DNA becomes part of the host's genetic makeup and is replicated.




# Target Plasmid

Vector : pET 28a




Enzyme site to use : Nde I, Xho I

Antibiotics : Kanamycin



# Template DNA information

Green fluorescent protein



# Template DNA information

Green fluorescent protein(GFP)

DNA sequence(750bp)

**ATG**GGCAGCAGCC**CATCATCATCATCAC**AGCAGCGGCCTGGTGCCGCGCGGCAGCCATATG**AAAGGAGAAGA**ACTTTT**CACTGGAGTTGTCCCAATTCTTGTTGAAT**  
**TAGATGGTGATGTTAATGGGCACGAATTTTCTGTCCGTGGAGAGGGTGAAGGTGATGCAACAATCGGAAA**ACTTACCCTTAAATTTATTTGCACTACTGGAAA**ACTACCT**  
**GTTCCATGGCCAACACTTGTCACTACTCTGACTTATGGTGTTC**AATGCTTTTCCCGTTATCCGGATCACATGAAACGGCATGACTTTTTCAAGAGTGCCATGCCCGAAGGT  
**TATGTACAGGAACGCACTATATCTTTCAAAGATGACGGGAAATACAAGACGCGTGCTGTAGTCAAGTTTGAAGGTGATACCCTTGTTAATCGTATCGAGTTAAAAGGTAC**  
**TGATTTTAAAGAAGATGGAAACATTCTCGGACACAACTCGAATACA**ACTTTAACTCACACGATGTATACATCACGGCAGACAAACAAGAAAATGGAATCAAAGCTGAA  
**TTCACTGTTCCGACAAACGTTGAAGATGGCTCCGTTCAACTAGCAGACCATTATCAACAAAATACTCCAATTGGCGATGGCCCTGTCTTTTACCAGACAACCATTACCT**  
**GTCGACACAACTGTTCTTTCGAAAGATCCCAACGAAAAGAGAGACCACATGGTCCTTCACGAGTACGTAACGCTGCTGGGATTACA** TGA

**Start codon**  
**Histidine tag**  
**GFP**  
**Stop codon**

## Primer

Forward primer (Nde I)

5' actg**CATATG**aaaggagaagaactttt**cactggagt** 3'

Reverse primer (Xho I)

5' actg**CTCGAG**tcatgtaatcccagcagc 3'

# Principles of PCR

## Organisms

### • A summary of DNA replication





Figure 16.16

Copyright © 2005 Pearson Education, Inc. publishing as Benjamin Cummings

### Essential components of DAN replication

1. Short complementary oligomer synthesized by **Primase**
2. **DNA polymerase**
3. **Helicase** for unwinding double stranded DNA
4. **Nucleotides**

## Laboratory



Copyright © 2004 Pearson Education, Inc., publishing as Benjamin Cummings.

### Essential components of PCR

1. **Primer**, which is artificially synthesized short oligomer
2. **Heat-stable DNA polymerase (pfu/Taq)**
3. **PCR machine** for the cycle of heating and cooling step
4. **Nucleotides** and **buffer** for polymerase

# Principles of PCR

## - Polymerase Chain Reaction(PCR)

Technology that amplifies DNA by repeatedly replicating specific DNA sequence. The DNA amplification time is short and a large amount of DNA can be obtained from very small amount of DNA.

## - Experiment term

Template DNA : Target DNA fragment to be amplified

Primer : Short oligonucleotide that binds to the part to be amplified

Pfu polymerase : Thermostable DNA polymerase. Low error rate

dNTP : A nucleotide that becomes a material for synthesizing a gene


Buffer : Contains various materials required for polymerase chain reaction

## - PCR step

Denaturation : Breaking hydrogen bonds in double helix. Double strand → Single strand

Annealing : Binding primers to DNA

Extension : Start synthesis of new DNA strands using dNTP



# Experimental process

- 1) Program as follows in PCR machine
- 2) Prepare reaction mixture as follows
- 3) Smooth pipetting, insert the PCR tube and run the program on the PCR machine.

Program setting

Cycle step	Temperature	Time	Cycle number
Initial denaturation	95°C	2min	1
Denaturation	95°C	30sec	25~35
Annealing	55°C	30sec	
Extension	72°C	30sec	
Final extension	72°C	5min	1
	4°C	hold	

Reaction mixture composition

Template	1 $\mu\text{l}$
Primer 1 (Forward)	1 $\mu\text{l}$
Primer 2 (Reverse)	1 $\mu\text{l}$
Pfu-DNA polymerase	1 $\mu\text{l}$
dNTP (10x)	2 $\mu\text{l}$
Buffer (10x)	2 $\mu\text{l}$
3 <sup>rd</sup> distilled water	12 $\mu\text{l}$
<b>Total volume</b>	<b>20 <math>\mu\text{l}</math></b>

\*\* PCR is very sensitive at a extremely small quantity of DNA. So, all reagents should be stored as aliquots, and tip, tube and water for PCR should be sterilized.

# Pipette Usage

- I. Use within a specific capacity range for each pipette
- II. Do not apply more force than necessary
- III. Works slowly when transferring solution
- IV. Do not place pipette upside down or lay down while containing solution
- V. Suction slowly from the top of the solution