


# **Exp2.**

# **Plasmid Prep and Restriction**


# Vector

- A vector is a DNA molecule used as a vehicle to transfer foreign genetic material into another cell.



- Origin of replication (ori)
  - : particular sequence in a genome at which replication is initiated
- Multiple cloning site (MCS)
  - : polylinker, is a short segment of DNA which contains many (up to ~20) restriction sites
- Selectable marker
  - : a gene introduced into a cell that confers a trait suitable for artificial selection
  - : often antibiotics resistance genes

# Vector : pET-21a




- Origin of replication (ori)

- Multiple cloning site (MCS)

- : Nde I , Xho I

- Selectable marker

- : Ampicillin (AmpR)



# Mini-prep

- A method of plasmid DNA extraction and purification from bacteria
- Buffer S1 (Glucose, Tris, EDTA, RNase A)
  - : open the cell wall, RNA degradation
- Buffer S2 (NaOH, SDS)
  - : open the cell wall, denatures chromosomal DNA and proteins, and releases plasmid DNA into the supernatant.
- Buffer G3 (Potassium acetate)
  - : renaturation, precipitates chromosomal DNA and protein
- Buffer PW (EtOH)
  - : removes salts and other cellular components
- Buffer EB
  - : plasmid elution

# Mini-prep protocol

- ① Pick up a single colony from fresh cultured LB agar plate (contains antibiotics) and inoculate the cell into the 5mL of fresh LB liquid media (contains antibiotics) at 37°C with shaking for 12~16hr.
  - Antibiotics : ampicillin (Amp)
- ② Transfer 1mL of bacteria culture to 1.5mL tube and centrifuge at 10,000rpm for 5 min at RT. Discard supernatant using pipette.
  - Labeling is needed. (Name or student number etc)
- ③ Resuspend pelleted bacteria cell thoroughly in 250uL of buffer S1 by pipetting until no clumps remain.
  - Buffer S1 contained RNase A must be stored at 4°C.
- ④ Add 250uL of buffer S2 and mix by inverting the tube 4 times (DO NOT VORTEX).
  - Buffer S1 is strong base, do not incubate for more than 3 min.
- ⑤ Add 350uL of buffer G3 and immediately mix by inverting the tube 4~6 times (DO NOT VORTEX).


# Mini-prep protocol

- ⑥ Centrifuge for 10 min.
- ⑦ Transfer carefully the supernatant to a spin column by pipetting. Centrifuge for 30 sec.
  - Labeling is needed. (Name or student number etc)
- ⑧ Remove the spin column, discard the pass-through, and re-insert the spin column to the collection tube.
- ⑨ Apply 700uL of buffer PW and centrifuge for 30 sec.
- ⑩ Remove the spin column, discard the pass-through, and re-insert the spin column to the collection tube.
- ⑪ Centrifuge for an additional 1 min to remove the residual wash buffer.
- ⑫ Transfer the spin column to a new 1.5mL microcentrifuge tube.
  - Labeling is needed. (Name or student number etc)
- ⑬ Add 50uL of buffer EB, let stand for 1 min, and centrifuge for 1 min.


# Restriction Enzyme (RE)

- Sequence-specific endonuclease
- Enzyme that cuts DNA at recognition site sequence (Palindrome)

- Sticky end (cohesive end)



- Blunt end



- 1 units
  - The amount of enzyme that can digest 1ug of DNA (50uL reaction volume) in 1 hr with optimal buffer condition

# RE digestion protocol

- ① Make follow reaction mixture.
  - RE should be on ice.
  - Well mixing by tapping and spin-down.

	Vector (mini-prep)	PCR product
DNA	16 uL	16 uL
Restriction buffer IV	2 uL	2 uL
Xho I (20U/uL)	1 uL	1 uL
Nde I (20U/uL)	1 uL	1 uL
Total vol.	20 uL	20 uL

- ② Incubate at 37°C for 2 hr.

