


**Exp6.**  
**DNA Purification**  
**(mini-prep, RE digestion, gel extraction)**


# Vector

- A vector is a DNA molecule used as a vehicle to transfer foreign genetic material into another cell.




- Origin of replication (ori)
  - : particular sequence in a genome at which replication is initiated
- Multiple cloning site (MCS)
  - : polylinker, is a short segment of DNA which contains many (up to ~20) restriction sites
- Selectable marker
  - : a gene introduced into a cell that confers a trait suitable for artificial selection
  - : often antibiotics resistance genes

# Vector : pET-21a



- Origin of replication (ori)
- Multiple cloning site (MCS)  
: Nde I , Xho I
- Selectable marker  
: Ampicillin (AmpR)



Ampicillin

# Mini-prep

- A method of plasmid DNA extraction and purification from bacteria
- Buffer S1 (containing Glucose, Tris, EDTA, RNase A) – 4°C storage (in fridge)
  - : open the cell wall, RNA degradation
- Buffer S2 (containing NaOH, SDS)
  - : open the cell wall, denatures chromosomal DNA and proteins, and releases plasmid DNA into the supernatant.
- Buffer G3 (containing Potassium acetate)
  - : renaturation, precipitates chromosomal DNA and protein
- Buffer PW (containing EtOH)
  - : removes salts and other cellular components
- Buffer EB
  - : plasmid elution

# Mini-prep protocol

- ① Pick up a single colony from fresh cultured LB agar plate (contains antibiotics) and inoculate the cell into the 5 mL of fresh LB liquid media (contains antibiotics) at 37 °C with shaking for 12~16hr.  
- Antibiotics : Ampicillin (Amp)
- ② Transfer 1mL of bacteria culture to 1.5 mL tube and centrifuge at 15,000 rpm for 5 min at RT.  
Discard supernatant using pipette.  
- Labeling is needed. (Name or student number etc)
- ③ Resuspend pelleted bacteria cell thoroughly in 250 uL of buffer S1 by pipetting until no clumps remain.  
- Buffer S1 contained RNase A must be stored at 4 °C.
- ④ Add 250 uL of buffer S2 and mix by inverting the tube 4 times (DO NOT VORTEX).  
- Buffer S1 is strong base, do not incubate for more than 3 min.
- ⑤ Add 350 uL of buffer G3 and immediately mix by inverting the tube 4~6 times (DO NOT VORTEX).


# Mini-prep protocol

- ⑥ Centrifuge for 10 min at 15,000 rpm.
- ⑦ Transfer carefully the supernatant to a spin column by pipetting. Centrifuge for 30 sec.  
- Labeling is needed. (Name or student number etc)
- ⑧ Remove the spin column, discard the pass-through, and re-insert the spin column to the collection tube.
- ⑨ Apply 700uL of buffer PW and centrifuge for 30 sec at 15,000 rpm.
- ⑩ Remove the spin column, discard the pass-through, and re-insert the spin column to the collection tube.
- ⑪ Centrifuge for an additional 1 min to remove the residual wash buffer.
- ⑫ Transfer the spin column to a new 1.5mL microcentrifuge tube.  
- Labeling is needed. (Name or student number etc)
- ⑬ Add 30 uL of buffer EB to the center of the membrane in the spin column.
- ⑭ Let stand for 1 min, and centrifuge for 1 min at 15,000 rpm.


# Restriction Enzyme (RE)

- Sequence-specific endonuclease
- Enzyme that cuts DNA at recognition site sequence (Palindrome)

- Sticky end (cohesive end)



- Blunt end




- 1 units
  - The amount of enzyme that can digest 1ug of DNA (50uL reaction volume) in 1 hr with optimal buffer condition

# RE digestion protocol

- ① Make follow reaction mixture.
  - RE should be on ice or cold rack.
  - Well mixing by tapping and spin-down.
- ② Incubate at 37°C for 1 hr.
  - Don't forget a labeling  
(your initial or student # etc)

	Vector (mini-prep)	PCR product
DNA	16 uL	16 uL
Restriction buffer IV	2 uL	2 uL
Xho I (20U/uL)	1 uL	1 uL
Nde I (20U/uL)	1 uL	1 uL
Total vol.	20 uL	20 uL





# Running an agarose gel

- ① Dissolve 1.5 g EB agarose powder to 100 mL 1 X TAE buffer.
  - Placing plastic wrap over the top of the flask to protect evaporation of buffer.
  - Microwave for 3 min until the agarose powder is completely dissolved.
- ② Let agarose solution cool down by stirring to about 50 °C.
- ③ Add 5 µL of SYBR safe.
  - SYBR safe allows you to visualize the DNA under blue light.
- ④ Pour the agarose into a gel tray with the comb, and let it until it has completely solidified.
- ⑤ Place the gel into the gel tank and fill the tank with 1 X TAE.
- ⑥ Mix the sample (RE digested vector and PCR product) with 2X loading dye.
  - Role of loading dye
    - 1) a visible dye - helps with gel loading & allows you to gauge how far the DNA has migrated;
    - 2) it contains a high percentage of glycerol → increases the density of your DNA sample causing it settle to the bottom of the gel well, instead of diffusing in the buffer.
- ⑦ Carefully load your samples (20 µL) and DNA marker (4 µL) into the wells of the gel.
- ⑧ Running gel for 130 V, 20 min.

# Gel extraction (Gel purification)

- A method to extract pure DNA fragment from agarose gel.
- Buffer GB (containing chaotropic salt)
  - : Enhance binding of DNA onto a membrane.
- Buffer NW (containing EtOH)
  - : removes salts and other cellular components
- Buffer EB
  - : DNA elution

# Gel extraction (Gel purification) protocol

- ① Visualize your DNA fragment using blue light.
  - Using the DNA ladder as a guide of the size of the DNA fragments
- ② Slice the desired DNA fragment from the gel using razor, and place it in a EP tube.
- ③ Add 400 uL of Buffer GB to the gel.
- ④ Incubate at 50 °C until the agarose gel is completely melted.
  - To help dissolve gel, mix by inverting the tube every 2-3 minutes.
- ⑤ Add 150 uL of isopropanol to the tube and mix.
  - To increase the yield.
- ⑥ Transfer the mixture to a spin column.
- ⑦ Centrifuge for 1 min at 6000 rpm. Discard the flow through and reinsert the column into the tube.
- ⑧ Add 750 uL of Buffer NW to the spin column. Let it stand for 5 min.
- ⑨ Centrifuge for 30 sec at 6000 rpm. Discard the flow through and reinsert the column into the tube.

# Gel extraction (Gel purification) protocol

- ⑩ Centrifuge for 1 min at 15000 rpm. Transfer the spin column to a new EP tube.
  - Please label your name at the side of the tube.
- ⑪ Apply 30  $\mu$ L of Buffer EB to the center of the membrane in the spin column.
- ⑫ Let it stand for 1 min, and centrifuge for 1 min at 15000 rpm.